(EN) Instructions for Use
(DE) Gebrauchsanweisung
(F) Mode d'emploi


(EN) **SedimSure 6% HES 600 in normal saline**
(DE) **SedimSure 6 % HES 600 in isotoner Kochsalzlösung**
(F) **SedimSure 6 % HES 600 en solution physiologique**

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(EN)

Introduction:

SedimSure is employed in the quantification of cellular components by sedimentation.

SedimSure is a ready-to use sedimentation aid solution for the addition to an anticoagulated harvest of hematopoietic stem cells for the purpose of volume reduction or red cell depletion by means of centrifugation according to the method described by Rubinstein et al (1) prior to cryopreservation. SedimSure is an aqueous solution of 6% HES 600 (hetastarch) and 0.9% sodium chloride (NaCl).

HES (hydroxy ethyl starch) enhances the erythrocyte sedimentation rate and in combination with adequate centrifugation enables an improved depletion of red blood cells from the blood fraction which contains the stem cell rich buffy coat.

SedimSure is a clear, pale yellow to amber solution. Exposure to prolonged adverse storage conditions may result in a change to a turbid deep brown or to the formation of a crystalline precipitate. Do not use the solution if these conditions are evident.

The container-solution system is a closed system and is not dependent upon entry of air during withdrawal of the solution. The container is overwrapped to provide an additional moisture barrier. The closure system has two ports for withdrawal by use of a syringe or a spike.

Exposure of SedimSure to heat should be minimized. Avoid excessive heat. Protect from freezing. It is recommended that the product should be stored at room temperature (25 °C); however, brief exposure up to 40 °C does not adversely affect the product.

The primary packaging container is sealed in an overwrap bag. This overwrap bag may only be removed ultimately prior to application of SedimSure. SedimSure is indicated for single use.

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Application of SedimSure

Use only if solution is clear and container and seals are intact.

Visual Checking

1. Do not remove the plastic container from its overwrap until immediately before use.
2. While the overwrap is intact, identify the solution SedimSure, the lot number, and the expiration date.
3. Check that the solution is clear.
4. Inspect the intact units for signs of obvious damage. If present, the unit should not be used.

Removal of Overwrap

To open overwrap, tear at any notch located at either end of unit. After removing the overwrap, check for any leakage by squeezing container firmly. If any leaks are found, discard unit as sterility may be impaired.

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Withdrawal of solution from container

1. Use aseptic techniques to avoid contamination of SedimSure and the target blood product as appropriate.
2. If possible use the hanger of the bag to hang it onto an appropriate holding device to position the bag vertically.
3. The SedimSure bag disposes of two access ports according to ISO 15747. The access ports are designated for puncture with a spike according to ISO 8536-4 or a cannula of 0.6 mm or 23 gauge outer diameter. After removal of the spike or cannula from the port the SedimSure bag does not form a sterile barrier anymore.
4. Twist off the plug from either port.
5. Insert the syringe needle respectively the spike via the port. If using a spike insert it with a twisting motion until the spike is firmly seated.
6. Withdraw solution aseptically and volumetrically and apply as described below.
7. Discard unit after single use.

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

According to the protocol of Rubinstein et al (1) 6% HES solution in normal saline is added to an anticoagulated cord blood harvest at a ratio of 1:5. The calculated amount required of SedimSure is withdrawn volumetrically, for example by means of a syringe or a spike from the SedimSure bag via one of the ports. The required amount is added to the anticoagulated target cell suspension. After addition the target solution is mixed gently to achieve a homogeneous distribution of the sedimentation aid throughout the cell suspension. The addition should be carried out under conditions beneficial for the cell suspension. The exact procedure of addition of SedimSure to the anticoagulated suspension of hematopoietic stem cells is to be validated by the user. After having obtained the volume reduced cord blood unit, a cryoprotectant is added. Then the cryopreservation process can start.

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Immediately after the completed thawing process SedimSure must be washed out of the stem cell suspension. This can be performed simultaneously in one procedure with the wash out of the cryoprotective solution. The washing process is executed in several washing steps consisting of centrifugation, withdrawal of liquid supernatant and resuspension of the cells with an appropriate wash solution. During the washing process until the quantitative elimination of SedimSure from the stem cell suspension, the suspension has to be kept cool at 2 °C. Consequently the wash out process has to be performed by means of a refrigerated centrifuge. The validation of the elimination process is the responsibility of the user. For performing the wash out process within a closed system several methods are described in scientific literature (2 - 8).

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Sources/Quellen/Bibliographie

- 1) Rubinstein, P.; Dobrila, L.; Rosenfield, R.E.; Adamson, J.W.; Migliaccio, G.; Migliaccio, A.R.; Taylor, P.E.; Stevens, C.E. : "Processing and cryopreservation of placental / umbilical cord blood for unrelated bone marrow reconstitution". Proc. Natl. Acad. Sci. USA 92, 10119-10122 (1995)
- 2) Calmels B, Houzé P, Hengesse JC, Ducrot T, Malenfant C, Chabannon C.: "Preclinical evaluation of an automated closed fluid management device: Cytomate, for washing out DMSO from hematopoietic stem cell grafts after thawing.", Bone Marrow Transplant. 2003 May;31(9):823-8.
- 3) Rodriguez L, Velasco B, Garcia J, Martín-Henao GA.: "Evaluation of an automated cell processing device to reduce the dimethyl sulfoxide from hematopoietic grafts after thawing.", Transfusion. 2005 Aug;45(8):1391-7.
- 4) Foïs E, Desmartin M, Benhamida S, Xavier F, Vanneau X, Rea D, Femand JP, Arnulf B, Mounier N, Ertault M, Lotz JP, Galicier L, Raffoux E, Benbunan M, Marolleau JP, Larghero J.: "Recovery, viability and clinical toxicity of thawed and washed haematopoietic progenitor cells: analysis of 952 autologous peripheral blood stem cell transplantations.", Bone Marrow Transplant. 2007 Nov;40(9):831-5. Epub 2007 Aug 27.
- 5) Laroche V, McKenna DH, Moroff G, Schierman T, Kadidlo D, McCullough J.: "Cell loss and recovery in umbilical cord blood processing: a comparison of postwash and postwash samples.", Transfusion. 2005 Dec;45(12):1909-16.
- 6) Lemarie C, Calmels B, Malenfant C, Arneodo V, Blaise D, Viret F, Bouabdallah R, Ladaïque P, Viens P, Chabannon C.: "Clinical experience with the delivery of thawed and washed autologous blood cells, with an automated closed fluid management device: Cytomate.", Transfusion. 2005 May;45(5):737-42.
- 7) Nagamura-Inoue T, Shioya M, Sugo M, Cui Y, Takahashi A, Tomita S, Zheng Y, Takada K, Kodo H, Asano S, Takahashi TA.: "Wash-out of DMSO does not improve the speed of engraftment of cord blood transplantation: follow-up of 46 adult patients with units shipped from a single cord blood bank.", Transfusion. 2003 Sep;43(9):1285-95.
- 8) Perotti CG, Del Fante C, Viarengo G, Papa P, Rocchi L, Bergamaschi P, Bellotti L, Marchesi A, Salvaneschi L.: "A new automated cell washer device for thawed cord blood units.", Transfusion. 2004 Jun;44(6):900-6.